

U-WAND / RECENT CONSUMER RESEARCH

Severe Acne Test

Two male volunteers with severe Acne were recruited for a pre-test where they followed the same regimen as the full trial volunteers for 2 weeks. Following is one of the before and after photo comparisons

PRE-TEST PROGRAM



Severe Acne Condition Before

Cleared of Acne after 2 Weeks

Unretouched Photos after Severe Acne test

Clinical Trial Demonstration & Results - 30 Days

To monitor the clinical trial, TSI used a company that is the industry leader in providing the imaging system and support necessary for high definition photographic documentation in the medical and skincare fields. According to their analysis, the average reduction in facial lesions from the use of the U-Wand Anti-Acne Treatment System was 70.3% over all the volunteers in 4-weeks. The highest reduction score achieved by this Company prior to this clinical trial was 20% reduction in lesion count.

